

R-AID Trolley PRO

Multipurpose backpack with trolley
Dimensional and structural change



Emergency trolley extremely versatile and robust, usable also as a backpack. 4 coaxial wheels ensure high load resistance while dragging the trolley on the ground.

Changes

Carrying handle moved

The handle, previously on the side of the backpack, is now on the front

Increased dimensions

The larger size, especially height, has allowed to increase the overall volume of 8 liters

Detachable backpack removed

There are no longer external removable parts

Integration of internal partitions

The main compartment contains dividers that allow to split the compartment from 2 up to 9 parts

Changes active from batch number 1501031

Accessories

4 R-Series internal bags with strap fixing

Technical data

Dimensions	620x350x410 ± 10 mm
PVC bag dimensions	360x160 mm
Weight	2,88 ± 0,2 kg
Volume ¹	48 ± 0,5 l
Shoulder strap length	From 70 to 130 cm
Trolley handle extension	43 cm
Materials	Nylon, PVC, PU
Loading capacity	30 kg

¹ Main compartment volume 32 l

Caratteristiche specifiche

Excellent portability thanks to the trolley system, 2 handles and shoulder strap

1 main padded compartment, equipped with straps to fix a bottle (max 2l) and removable internal dividers that allow to split the compartment up to 9 parts

1 transparent pocket with zip closure on the internal side of the closure flap of the main compartment

2 Transparent pockets with zip closure of which:

- 1 equipped with internal pocket with zip closure

- 1 equipped with *transparent detachable bag* with zip closure, fixed by means of straps

1 Front transparent pocket with strap closure

1 Front pocket with zip closure

1 top pocket with zip closure and carrying handle

Padded shoulder straps and lumbar belt with quick release buckle improve load distribution

Reflective inserts for increased visibility

Reinforced bottom

R-AID TROLLEY PRO - RED/BLACK CORDURA TROLLEY
CB04180B

Rev.0 (21/12/2015)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it