


Structural modification
RESPONSE SERIES
 Multi-use emergency bag



RESPONSE



FIRST RESPONSE



PRO RESPONSE

A revision has been made to the number of acetylic rings for shoulder strap fixture and the number of shoulder straps supplied.

The Response Series bags: Response, First Response and Pro Response will be supplied with just one shoulder strap.

The rings for its fixture will be positioned diagonally on each side of the bag so that the bag is perfectly balanced even when loaded.

TECHNICAL DATA

Dimensions (mm):	465 x 270 x h270
Capacity (L):	34
Weight (kg):	1,8
Material:	Oxford 840 D
Carry strap:	Yes
Number of handles:	1
Rigid double bottom:	Yes
Bottom reinforcement:	Yes
Reflective insert:	Yes
Cylinder fixation:	Yes
Internal material:	Heat moulded plastic

Personalization methods

Printing space (mm):	100 x 50
Part to be personalized:	Top flap

Response orange	CB07010A
Response red	CB07011A
Response blue	CB07012A
First Response	CB07014A
Pro Response	CB07021A

The revision described is valid from 15/09/2011