

Structural revision and upgrade

Fixo Splint

Rigid splints



QC70050A



QC70020A



QC70010A



Both the shape and the production material of all three sizes have been revised and the product is now notably improved when compared with the original.

The new shape offers improved strength and the welded edges avoid the penetration of any liquid inside and make the clearing process more efficient so that cross infections risks are greatly reduced.

TECHNICAL DATA

Length size. XS	310 mm
Length size S	510 mm
Length size. L	740 mm
External material	Nylon coated PVC
Internal material	Closure
Temperature of use	From 0 to +40 °C
Temperature for storage	From -15 to +60 °C

QC70010A	Fixo Splint size XS
QC70020A	Fixo Splint size S
QC70050A	Fixo Splint size L
QC70090A	Fixo Splint Kit Splint sizes S e L
QC70091A	Fixo Splint Kit Splint sizes XS, S e L

The revision described is valid from 31/10/2013