

Product revision

RES-Q-SPLINT E RES-Q-SPLINT PLUS

Transport bag for vacuum splints



The special transport bag for the Res-Q-Splint e Res-Q-Splint Plus vacuum splints will undergo some revisions in the next production.

The bag will no longer be sewn and to increase the product strength and durability it will be welded ultrasonically. The colour will no longer be orange. It will be yellow for increased visibility.

A new logo will make the identification of the emergency immobilising devices inside faster

QM22540A – Transport bag for Res-Q-Splint e Res-Q-Splint Plus vacuum splint

TECHNICAL DATA

Dimensions	550 x 180 x h420 mm
Weight	600 g
Material	Spentex®

The revision described is valid from 19/03/2014