

Revised products' configuration

ECO B-LIFE

Rubber face masks



When stock terminates the rubber face masks Eco B-life will be removed from production. In configurations of manual resuscitators which included the presence of this model of face masks and in resuscitation kits inside where they were until now present, the Eco B-life will be substituted with the following codes:

- RM20400B - B-life Silicone face mask size 0
- RM20401B - B-life Silicone face mask size 1
- RM20402B - B-life Silicone face mask size 2
- RM20403B - B-life Silicone face mask size 3
- RM20404B - B-life Silicone face mask size 4
- RM20405B - B-life Silicone face mask size 5
- RM20410B - B-life Silicone face masks kit six sizes

The B-life masks are characterized by a silicone cushion, which improves the fit of the mask to the patient's face. The correspondences between the sizes remain unchanged.

DATI TECNICI

Material of shell	polycarbonate
Material of cushion	silicone

Sterilization with autoclave	at 134°C
------------------------------	----------

ENVIRONMENTAL CONDITIONS

During functioning	from -18 to +50 °C
Durante storage	from -30 to +60 °C

Latex free

CE 0123

The revision described is valid from 15/05/2014