

Revised configuration

B-VAC

Manual aspirator for secretions



The B-Vac suction device will now be supplied complete with only the adult catheter.

The paediatric catheter will no longer be included in the set.

TECHNICAL DATA

Dimensions (mm)	80 x 70 x 35
Weight (kg)	0,29
Maximum vacuum (mmHg)	-550
Maximum capacity (L/min)	>20
Volume supplied container	350
PVC tube	Yes
Temperature of use (°C)	-20 +50
Temperature for storage (°C)	-10 +60
Humidity	95%
Color	Yellow

B-Vac RV71150B

The revision described is valid from 12/03/2013