

Structural modification
SPENCER SHELL
Universal basket stretcher



Thanks to a modification of the structure material the lifetime of the Spencer basket stretchers Shell has been doubled from 5 to 10 years. The polyethylene now being used has a new formulation that presents an increase of the mechanical performances and of the rigidity that permit the raising of the Shore values. This thanks to new systems of polymerization of the raw material. The modification of the material confers the device a prolonged lifetime while it completely maintains the Spencer security and quality requirements.

TECHNICAL DATA

Spencer Shell yellow	ST04012A
Spencer Shell orange	ST04000A
Spencer Shell red	ST04017A
Spencer Shell military	ST04001A



The revision described is valid from 06/06/2012