

Functional modification

## R-MAX

10 G Fixing system for stretchers



The modification concerns the release lever of the hook, which has been modified to make the activation easier and more ergonomic.

Compared to the previous version, now the operator has no longer to pull the release lever, but simply turn it toward the outer side of the hook, moreover the new internal mechanism reduces the force required to activate it.

The system is equipped with a new ergonomic grip suitable for all operators.

### TECHNICAL DATAS

Length	1800 mm
Width	445 mm
Height	75 mm
Weight	9 kg

ST42702C

R-MAX



The revision described is valid from 19/04/2012