

Technical and functional upgrade

EXTRATANK

Recovery chest



The excellent quality of the Spencer Extratank recovery chest has now undergone further improvements with structural alterations in specific areas. The new aluminium frame with steel reinforcements makes the device extremely rigid and lightweight. Comfort for the operator is increased thanks to the side transport handles and the telescopic levers with rubber covering.

TECHNICAL DATA

Upper shell	High density polyethylene
Lower shell	High density polyethylene
Frame	Steel/aluminium
Wheels Ø 200 [mm]	Rubber polyurethane
Telescopic handles	Aluminium
Folding handles	Steel
Closure	Elastic rope
Width [mm]	640
Length with closed handles [mm]	2240
Length with open handles [mm]	2770
Height [mm]	510
Volume [L]	400
Weight [kg]	25
Load capacity [kg]	220

Extratank grey ZC00201A
Extratank white ZC00204A

The revision described is valid from 03/10/2012